**Actim Combi**® Hemoglobin. Transferrin.

## INSTRUCTIONS FOR USE

- |                                |                                    |
|--------------------------------|------------------------------------|
| <b>CS</b> Návod k použití      | <b>IT</b> Istruzioni d'uso         |
| <b>DA</b> Brugsvejledning      | <b>NL</b> Gebruiksaanwijzing       |
| <b>DE</b> Gebrauchsanweisung   | <b>NO</b> Bruksanvisning           |
| <b>ES</b> Instrucciones de uso | <b>PT</b> Instruções de utilização |
| <b>FI</b> Käyttöohje           | <b>SV</b> Bruksanvisning           |
| <b>FR</b> Notice d'utilisation |                                    |

Medix  
Biochemica





[www.medixbiochemica.com](http://www.medixbiochemica.com)

## STRUCTURE OF DIPSTICK

## TEST PROCEDURE AND RESULTS

# INSTRUCTIONS FOR USE

Numbers ① - ⑧ refer to illustrations on inner cover.

## STRUCTURE OF DIPSTICK

① Dip area ② Result area ③ Test lines ④ Control line

## INTENDED USE

The Actim Combi<sup>®</sup> Hemoglobin Transferrin test is a visually interpreted, qualitative immunochromatographic dipstick test for detecting occult blood in feces. The Actim Combi Hemoglobin Transferrin test detects human hemoglobin and human transferrin in a feces suspension. The test is intended for professional use to help diagnose gastrointestinal bleeding.

## KIT COMPONENTS

The components of the Actim Combi Hemoglobin Transferrin test are:

- Actim Combi Hemoglobin Transferrin dipstick in a perforating device. The dipsticks [10] are packed in an aluminum container.
- Specimen Dilution Buffer [3 ml]. This phosphatebuffered solution contains bovine serum albumin [BSA], protease inhibitors and <0.1 % (w/v) sodium azide [NaN<sub>3</sub>]. One sampling stick is attached to the cap of each tube.

The kit 30322ETAC contains 2 aluminum tubes of dipsticks [30351ETAC] and 20 Specimen Dilution Buffer tubes [30371ETAC].

An instructions for use leaflet is included in all product formats. Toilet paper or a clean disposable container is needed for sampling (not provided with the test).

## STORAGE

Store the dipsticks and the Specimen Dilution Buffer tubes at +2...+30 °C. Stored unopened, each component can be used until the expiry date marked on the component. After opening the aluminum container first time, the dipsticks can be stored for 4 months at +2...+30 °C. However, the expiry date must not be exceeded.

## SPECIMEN COLLECTION

The specimen is feces suspension that is prepared as soon as possible after sample collection. The original feces sample can be stored for 3 days at +2...+8 °C.

Unscrew the cap of the Specimen Dilution Buffer tube. Collect the feces sample by using the sampling stick attached to the screw cap. There are two small slits in the head of each stick. Twist the stick randomly in several different places in the feces sample so that both slits become filled with feces. In case of liquid samples, dry the stick by wiping it with a piece of e.g. toilet paper before sticking it into the sample. The slits of the stick will be filled with the liquid feces.

Put the sampling stick into the Specimen Dilution Buffer tube by pushing the stick through the cone. Close the cap tightly and shake the tube to suspend the feces in the buffer. The specimen suspension is now ready for testing. The suspension can be stored for one week at +2...+25 °C (+2...+8 °C is recommended).

## TEST PROCEDURE AND INTERPRETATION OF THE RESULTS

1. If stored refrigerated, allow the aluminum container of dipsticks and the buffer tube to reach room temperature before opening. Remove the required number of perforating devices (each containing a dipstick) from the container just prior to use and close the container immediately. Do not touch the yellow dip area at the lower ends of the dipsticks. Do not touch the part of the perforating device below the stopper, as this area will come into contact with the sample. Identifying marks may be written on the upper red parts of the dipsticks.
2. The specimen should be at room temperature upon testing. Use the sharp tip of the perforating device to puncture the perforation area in the bottom of the buffer tube **then push the device into the tube until the stopper prevents it going further** ⑤.
3. Slowly invert the tube upside down once and immediately place the tube on the table in an upright position ⑥. **NOTE! It is important that the inverting movement is not too abrupt, a peaceful inversion is vital for the correct performance of the test. Also, the tube must not remain upside down.**
4. The result can be interpreted as positive as soon as one blue line or one red and or both blue and red lines, and one black line become visible in the result area. The negative result should be read at **10 minutes** after the liquid front has reached the lower part of the result area ⑦. **Do not pay attention to any lines appearing later than 10 minutes.**

#### 5. If, in addition to the black control line you see:

- **one blue line** in the **lower part** of the result area, the result is **hemoglobin positive**.
- **one red line** in the **middle part** of the result area, the result is **transferrin positive**.
- **one blue line and one red line** in the result area, the result is **hemoglobin and transferrin positive**.
- **no blue or red lines** in the result area, the result is **negative**.

If **the control line does not appear**, the test is **invalid** Ⓢ.

In case of a missing control line, please see the Notes section for troubleshooting.

#### DISPOSAL OF THE TEST DEVICE

The contact between the Specimen Dilution Buffer tube cap and the perforating device is so close that the test device can be disposed of directly after the result has been read. Do not remove the perforating device from the tube.

#### LIMITATIONS OF THE TEST

The test is intended for *in vitro* diagnostic use only.

#### NOTES

- The test results are qualitative. No quantitative interpretation should be made based on the test results.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick if you notice a blue coloring in the result area before testing.
- Take care that any part of the test device that will be in contact with the sample tube do not become contaminated by external blood, e.g. from a cut finger. Do not touch the dip area of the dipstick or the area of the perforating device below the stopper.
- When inverting the sample tube to allow the dipstick to absorb sample, do not invert the tube too fast, or for too long. The test will not work properly if the amount absorbed is too small or too large.
- If the control line does not appear, the inversion of the sample tube may have

been performed too abruptly. Take a new dipstick and repeat the inversion in a controlled way.

- Improper sampling may lead to false result.
- Occult blood is not uniformly distributed in the feces sample. Therefore, it is important to dip the sampling stick randomly in several places in the feces sample.
- Appearance of a control line confirms correct performance of the test. If a control line does not appear the test is invalid, and should be repeated using another dipstick.
- If the test result cannot be interpreted clearly (e.g. if the lines are blotched or uneven) it is recommended that the test be repeated with another dipstick.
- At ten minutes the appearance of any faint-to-dark blue and/or red test line along with a control line indicates a positive result. However, do not pay attention to any lines appearing after 10 minutes.
- If only the control line is visible, the test result should be interpreted as negative only after 10 minutes have elapsed.
- If the color of the test line is not blue (hemoglobin) or red (transferrin), the result should be interpreted as invalid and the test should be repeated with a new feces sample.
- Patients with the following conditions should not be considered for testing, as these conditions may interfere with test results: bleeding hemorrhoids, constipation bleeding, menstrual bleeding and urinary bleeding. However, these patients may be considered for testing after such bleeding ceases.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

## PRINCIPLE OF THE TEST

Feces do not normally contain measurable amounts of blood. The presence of hemoglobin and transferrin in feces implies that there may be bleeding in the gastrointestinal tract.

The test is based on immunochromatography. It involves antibodies to human hemoglobin and to human transferrin. Antibodies are bound to blue (hemoglobin) or red (transferrin) latex particles (the detecting label). Antibodies against both types are immobilized on a carrier membrane to catch the complex of antigen and latex-labeled antibody and indicate a positive result.

When the feces suspension comes into contact with the dip area of the dipstick, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains hemoglobin and /or transferrin, they bind to the antibodies attached to the latex particles. The particles are carried by the liquid flow and, if hemoglobin and / or transferrin are bound to them, they bind to the corresponding catching antibodies. A blue line will appear in the result area if the concentration of hemoglobin in the sample exceeds the detection limit of the test. A red line will appear in the result area if the concentration of transferrin in the sample exceeds the detection limit of the test. A black line (control line) confirms correct performance of the test.

## PERFORMANCE OF THE TEST

### Analytical Sensitivity

The analytical sensitivity, detection limit, of the Actim Combi Hemoglobin Transferrin test was evaluated by two operators using samples with different concentrations ranging of human hemoglobin from 10 – 250 000 µg/l and human transferrin from 1 – 10 000 µg/l on three different lots of Actim Combi Hemoglobin Transferrin test.

The detection limit of the Actim Combi Hemoglobin Transferrin test is approximately 50 µg/l of human hemoglobin in the sample and the results remain positive at least to 250 000 µg/l. A feces suspension with 50 µg of hemoglobin/l corresponds to 25 µg of hemoglobin/g of feces.

The detection limit of the Actim Combi Hemoglobin Transferrin test is approximately 5 µg/l of human transferrin in the sample and the results remain positive at least to 10 000 µg/l. A feces suspension with 5 µg/l of transferrin corresponds to 1.5 µg of transferrin /g of feces. Transferrin 2.5 µg/l can be positive or negative which corresponds to 0.75 µg/g of feces.

### Repeatability and reproducibility

Intra- (repeatability) and inter- (reproducibility) assay precision were evaluated using samples containing 0 – 200 µg/l of human hemoglobin and transferrin in samples. For repeatability the samples were tested by two operators on the same day using 10 replicates per sample on three different lots of Actim Combi Hemoglobin Transferrin, the total amount being 140 tests/Lot (420/3 lots) in one day. For reproducibility the samples were tested on the five consecutive workdays using three replicates per sample on three different lots of Actim Combi Hemoglobin Transferrin, the total amount being 210 tests/Lot (630/3 lots) in five days. Repeatable and reproducible results were obtained.

### Analytical specificity, Cross-reactivity

The Actim Combi Hemoglobin Transferrin test was tested with following cross-reactivity's: Cross-reactivity of human transferrin in hemoglobin-line, cross-reactivity of human hemoglobin in transferrin-line, cross-reactivity of Bovine hemoglobin, Porcine hemoglobin, Bovine transferrin, Porcine transferrin and Human Lactoferrin. No cross-reactivity was observed. The Actim Combi Hemoglobin Transferrin test is specific for the human hemoglobin and transferrin.

### Method comparison

The function of Actim Combi Hemoglobin Transferrin test was compared to the immunochromatographic test of CerTest FOB-Transferrin. The agreement with the two methods was 97.1 % (hemoglobin) and 93.3 % (transferrin). The results are shown in FIG 1 on the inner back cover.

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FIG 1.


**COMPARISON RESULTS OF THE HEMOGLOBIN LINE:**

Actim Combi® Hemoglobin Transferrin	Reference: CerTest FOB-Transferrin		
	POS	NEG	TOTAL
POS	26	3	29
NEG	0	76	76
<b>OVERALL AGREEMENT 97.1%</b>			n=105

**COMPARISON RESULTS OF THE TRANSFERRIN LINE**

Actim Combi® Hemoglobin Transferrin	Reference: CerTest FOB-Transferrin		
	POS	NEG	TOTAL
POS	35	3	38
NEG	4	63	67
<b>OVERALL AGREEMENT 93.3%</b>			n=105

**EXPLANATION OF SYMBOLS** • Vysvětlení symbolů • Symbolforklaring • Interpretation der Symbole • Significado de los símbolos • Symbolien selitykset • Explication des symboles • Spiegazione dei simboli • Interpretatie van de symbolen • Symbolforklaring • Explicação dos símbolos • Teckenförklaringar

YYYY-MM-DD

Use By  
Použitelné do  
Anvendes inden  
Verwendbar bis  
Fecha de caducidad  
Viimeinen käyttöpäivä  
Utiliser jusque  
Utilizzare entro  
Houdbaar tot  
Brukes før  
Válido até  
Använd före

Temperature limitation  
Teplotní rozmezí od do  
Temperaturbegrensning  
Temperaturbereich  
Limite de temperatura  
Lämpötilarajoitus  
Limites de température  
Limiti di temperatura  
Temperatuurlimiet  
Temperaturbegrensning  
Limite da temperatura  
Temperaturbegrensning

Manufacturer  
Výrobce  
Producent  
Hersteller  
Fabricante  
Valmistaja  
Fabricant  
Fabbriicante  
Fabrikant  
Produsent  
Fabricante  
Tillverkare

**LOT**

Batch code  
Číslo šarže  
Batchnumber  
Chargenbezeichnung  
Codigo de lote  
Erä  
Code du lot  
Codice del lotto  
Lot nummer  
Lotnummer  
Número de lote  
Lot nummer

**REF**

Catalogue number  
Katalogové číslo  
Katalognummer  
Bestellnummer  
Número de catálogo  
Tuotenumero  
Référéncia du catalogue  
Numero di catalogo  
Catalogusnummer  
Katalognummer  
Código do produto  
Katalognummer

Single use  
Jednorázové použití  
Engangsbrug  
Nicht wiederverwendbar  
No reutilizar  
Kertakäyttöinen  
Ne pas réutiliser  
Non riutilizzare  
Niet opnieuw gebruiken  
Ikke til gjenbruk  
De uso único  
Återanvänds ej

Contains sufficient for <n> tests  
Lze použít pro <n> testů  
Indeholder tilstrækkeligt til "n" test  
Ausreichend für "n" Tests  
Contenido suficiente para "n" ensayos  
Sisältö on riittävä "n" testiä varten  
Contenu suffisant pour "n" tests  
Contenuto sufficiente per "n" saggi  
Inhoud voldoende voor "n" testen  
Inneholder nok til "n" tester  
Conteúdo suficiente para <n> exames  
Räcker till "n" antal tester

**DILSPE**

Specimen dilution buffer  
Pufř k ředění vzorku  
Fortyndingsbuffer  
Probenverdünnungspuffer  
Tampón de dilución de Muestra  
Näytteenlaimennuspuskuri  
Tampon de dilution  
Tampone di diluzione  
Buffer om de monsters te verdunnen  
Fortynningsbuffer  
Tampão para diluição da amostra  
Provutspädningslösning

**STICK**

Dipstick  
Testovací proužek  
Teststrimmel  
Teststreifen  
Tira reactiva  
Testitikk  
Bandelette  
Dipstick  
Dipstick  
Teststrimmel  
Tira reactiva  
Teststicka

Consult instructions for use  
Viz návod k použití  
Se brugsanvisningen  
Gebrauchsanweisung zurate ziehen  
Consultar las instrucciones de uso  
Katso käyttöohjeet  
Consulter la notice d'utilisation  
Consultare le istruzioni d'uso  
Raadpleeg de gebruiksinstructies  
Se bruksanvisningen  
Consultar as instruções de utilização  
Läs i bruksanvisningen

**IVD**

*In Vitro* Diagnostic Medical Device  
*In Vitro* diagnostický zdravotnický prostředek  
Medicinsk udstyr til *in vitro*-diagnostik  
*In vitro* Diagnostikum  
Producto sanitario para diagnóstico *in vitro*  
*In vitro* -diagnostiikkaan tarkoitettu  
lääkinnällinen laite  
Dispositif médical de diagnostic *in vitro*  
Dispositivo medico-diagnostico *in vitro*  
Medisch hulpmiddel voor *in-vitro*diagnostiek  
Medisinsk utstyr for *in vitro* diagnostikk  
Dispositivos médicos de diagnóstico *in vitro*  
Medicintekniska produkter för *in vitro* diagnostik

**CONT NaN<sub>3</sub> <0.1%**

Contains <0.1 % NaN<sub>3</sub>  
Obsahuje <0,1 % NaN<sub>3</sub>  
Indeholder <0,1 % NaN<sub>3</sub>  
Enthält <0,1 % NaN<sub>3</sub>  
Contiene <0,1 % NaN<sub>3</sub>  
Sisältää <0,1 % NaN<sub>3</sub>  
Contient <0,1 % NaN<sub>3</sub>  
Contiene <0,1 % NaN<sub>3</sub>  
Bevat <0,1 % NaN<sub>3</sub>  
Inneholder <0,1 % NaN<sub>3</sub>  
Contém <0,1 % NaN<sub>3</sub>  
Innehåller <0,1 % NaN<sub>3</sub>

# Actim Combi® Hemoglobin. Transferrin.

## Medix Biochemica

Oy Medix Biochemica Ab  
Klovinpellontie 3, FI-02180 Espoo, Finland  
Tel. +358 9 547 680  
[www.medixbiochemica.com](http://www.medixbiochemica.com)  
[medix@medixbiochemica.com](mailto:medix@medixbiochemica.com)

© Medix Biochemica. All rights reserved.  
Actim is a trademark of Medix Biochemica.

OACE30322-1 03/2017